## User Instruction Manual

### Oxford® Journey

To avoid injury, read user manual prior to use. For alternative languages, contact your authorised service provider. Always refer to [www.oxfordhealthcare.co.uk](http://www.oxfordhealthcare.co.uk) for latest revision.

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

# Contents


- 1. The Oxford Journey Patient Lift .....3
- 2. Introduction: About Your Lift .....4
- 3. Assembly & Commissioning Instructions .....5
- 4. Sling Guide.....10
- 5. Lifting With Your Oxford Journey.....12
- 6. Smart Monitor & Handset .....14
- 7. Operating Instructions .....18
- 8. Safety Precautions .....21
- 9. Charging Instructions.....23
- 10. Maintenance Schedule & Daily Check List .....25
- 11. Technical Specifications .....27
- 12. Warranty .....30

# 1. The Oxford Journey Patient Lift

## Ready for use position

## Storage position

## 2. Introduction: About Your Lift

Each Journey lift is fully assembled, load tested and certified before being packed/shipped.

The packing consists of a strong, purpose built carton to ensure the safe arrival of the lift. A number of documents are supplied in a wallet, and packed with each lift and should be kept safely for future reference:

- TEST CERTIFICATE
- PACKING CHECK LIST
- USER MANUAL

The TEST CERTIFICATE is an important document and will be required for your insurance records. It is valid for six months but after it has expired the lift should be inspected and serviced as per the maintenance schedule.

Your authorised dealer can carry out servicing and periodic testing for those countries where this is a requirement (Europe). Please ensure your lift is included in their maintenance schedule. If you are at all unsure what your local market servicing requirements are, please check with your dealer, a local government agency and/or Oxford.

### WARNING

**Prior to each daily use of your lift the checks detailed in the daily checklist within this manual should be carried out.**

**The CE mark:** The Oxford Journey carries the CE mark and complies with the following EC directives:

- Medical Device Regulation 2017/745
- EMC Directive 2014/30/EU
- Low Voltage Directive 2014/35/EU

### Statement of Intended Use (Journey)

The intended use of this lifting device is for the safe lifting and transfer of an individual from one resting surface to another (such as a bed to a wheelchair). Oxford recommends that the transfer of a patient is fully risk assessed and conducted safely over a short distance only.

The Oxford Journey is suitable for patients in the SITTING position only who have a degree of weight-bearing ability but require assistance to stand. By removing the foot tray, the Oxford Journey patient lift can also be used as a walking aid following a careful risk assessment.

The Oxford Journey is designed to support and promote safe patient handling and transfer for both the patient and carer. It can be folded without the use of tools enabling it to be easily stored in environments where space is limited or to facilitate onward transportation in a vehicle with a suitably sized boot.

The Oxford Journey stand aid is designed to be used in conjunction with the Oxford range of slings. The examples of slings suitable for use with this device are listed as follows:

- Oxford Standing Sling
- Oxford Transport Sling

### WARNING

**Oxford Patient Lifts should only be operated by fully trained/competent care givers.**

It is important that users of this product are adequately trained in correct use and safe manual handling practice. Oxford strongly recommends a face-to-face patient risk assessment is conducted to determine the safety and well-being of the patient and carer(s) at all times, before each and every use. Oxford cannot be held responsible for incorrect or untrained use that results in injury to the patient and/or carer. If you require training on this or any other product from Oxford, please contact us for support.

## Expected Service Life

Oxford Professional lifts are designed and tested for a minimum service life of ten (10) years, subject to the use and maintenance procedures stated in this manual. Use, other than in accordance with these instructions, may compromise service life.

## Serious Incident Reporting

In the event of a serious incident taking place during use of this product, affecting the patient and/or care-giver, it must be reported to the product manufacturer or authorised distributor. Should the incident take place within the European Union (EU), it must also be reported to the local competent authority within the member state.

# 3. Assembly & Commissioning Instructions

## Carton Contents

Place the carton in a clear working area and open carefully. The carton contains:

- OXFORD JOURNEY LIFT
- WALLET CONTAINING DOCUMENTS
- HAND CONTROL
- BATTERY PACK
- CHARGING LEAD
- DESK TOP CHARGER/STAND

## WARNING




**The Oxford Journey is heavy and will need to be lifted with care. You may need assistance to lift the Oxford Journey from the carton.**

## Assembly Instructions


The Oxford Journey is delivered in its folded/storage position; the following assembly process demonstrates the 3 easy steps required to take the Oxford Journey from its storage position to its ready for use position. This process is carried out with no additional tools required.

1. Remove all the parts from the carton and place on the floor, taking care to protect the finish from damage.

Lift the Oxford Journey carefully from the carton. Lie the lift flat on the floor and ensure the rear castors are securely locked.


2. To unfold the Oxford Journey into its ready-to-use position, raise the mast assembly upright. This is achieved by pressing the mast release button located under the actuator (as shown). This will release the mast assembly and allow the mast to be raised. Using the top of the push handle, raise the mast to its upright position and locate the mast into the slot (1). Once the mast is firmly located, use the mast locking lever to secure the mast. This is achieved by pushing the mast locking lever down towards the mast surface (2). Finally, release braked castors at the rear and the Oxford Journey is ready for use (3).

### **⚠ CAUTION**

Avoid trapping fingers. Keep fingers away from the end of the mast when inserting into the mast slot. Tighten the mast-locking device, and unlock the braked castors at the rear.

3. Adjust the cow horn mechanism into one of three positions to suit either the patient task and/or environment. To adjust the cow horn, simply pull the trigger backwards towards the rear of the lift while holding the cow horn. Then select the desired position as indicated on the positioning decal (1) and release the trigger to lock. The lift is now ready for operation.


**NOTE:** See page 12 for cow horn position guidelines.


## Disassembly

To disassemble and store the Oxford Journey follow the four easy steps shown below; please note that no additional tools are required.


1. Lower the boom into its lowest possible position and adjust the cow horn to the storage/fold position (as shown).


2. Ensure the rear castors are securely locked (1). Raise the knee pad to its upper most position (2).


3. Unlock the mast-locking lever by pulling it outwards from the mast (1). Using the push handle lift the mast upwards and then slowly push downwards (2), carefully folding the mast assembly towards the legs until you hear a “Click” (3). The “Click” indicates the mast is now in its locked position.

4. When folded, lift the Oxford Journey upright using the legs and push handle as a guide.

**NOTE:** If you wish to separate the mast and boom completely from the base and legs, undo the mast-locking lever fully (remove pin assembly from mast) and lift the mast and boom completely free from the base. Be careful not to lose the mast lock lever and location spindle.


**⚠ CAUTION**


Care should be exercised when folding/unfolding the Oxford Journey as there is a possible danger of trapping the fingers etc. Please follow the instructions carefully and ask for assistance if you are unsure of the correct procedure.

**ALWAYS CHECK THE FOLLOWING BEFORE OPERATION:**


- The mast is fully locked into position (see assembly instructions, page 6).


- The legs of the lift open and close freely (Using foot pedals).


- The power pack is fitted to the lift and the latch holding the power pack in place is fully engaged. A "Click" should be heard when it is in place.

- The red emergency stop button, located on the control box, is in the RUN (out) position.

- The up and down buttons on the hand control lower and raise the boom correctly.

## 4. Sling Guide


The Oxford Journey is suitable for patients in the SITTING position. The slings suitable for this device are listed as follows:

- \* Oxford Transport sling (various sizes available)
- \* Oxford Standing sling (various sizes available)

### USER GUIDE FOR STANDAIID SLINGS

Both the Transport and Standing sling are easy fit slings, suitable for standing and supporting. They are not a general purpose sling. Patients must always be assessed for suitability by a qualified person.

#### Transport Sling


#### Fitting Instructions

Feed the sling down the back of the patient leaving the top of the commode aperture at the base of the spine. Check the roll part of the sling is square across the patient's back. Attach the support strap around the chest of the patient. This strap need not be tight. It is there to hold the sling in position while attaching the loops to the Standaid. It may be left off altogether if preferred. Ensure patient's arms are on the outside of the sling (refer to Fig. C, page 13).

Raise the patient's leg and feed the leg strap under and up between the legs. Ensure the sling is not twisted or creased under the thigh. Repeat this procedure for the other leg. Attach the loops on the roll parts of the sling to the white hooks at the end of the handgrips. Without crossing the loops on the leg straps, attach to left and right hooks on the lifting fork. You would normally use the middle of the three loops. The patient is now ready to be raised from the seat. Once clear of the seat, the patient will be more comfortable if lowered back in to the sitting position. Oxford recommends that slings be checked regularly and particularly before use for signs of fraying or damage. Do not use slings that are worn or damaged.

### Standing Sling

### Fitting Instructions

Place the sling down behind patients back between base of shoulder blades and bottom of rib cage. Attach the support straps around the front of the patient. This strap need not be tight. It is there to hold the sling in position whilst attaching the sling to the Standaid. It may be left off altogether if preferred. Ensure that the patient's arms are outside the sling (refer to Fig. A, page 13).

**NOTE:** For detailed fitting instructions, please refer to the user guide supplied with each sling.

### **⚠ WARNING**

Oxford recommends that slings be checked regularly and particularly before use for signs of fraying or damage. **DO NOT** use slings that are worn or damaged.

### **⚠ WARNING**

**OXFORD RECOMMENDS THE USE OF GENUINE OXFORD PARTS.** Oxford sling and lift products are designed to be compatible with one another. For country specific guidance on sling use and compatibility, please refer to the sling label or contact your local market distributor or Oxford.


### **⚠ WARNING**

Refer to maximum weight capacity of lift. Sling capacity is limited by the maximum capacity of the lift.

## 5. Lifting With Your Oxford Journey


The Oxford Journey has an adjustable cow horn. This allows the carer to position the cow horn in a comfortable position to suit the patient or environment. The cow horn has 3 positions. Position 1 is more suitable for a large/taller patient, position 2 is more suitable for a medium/average size patient and position 3 is more suitable for a shorter/pediatric patient or when standing a patient from a low chair. (Please note, these are guidelines only, it may be necessary to adjust according to patient comfort). To adjust the cow horn, simply pull the trigger backwards towards the rear of the lift whilst holding the cow horn. Then select the desired position as indicated by the positioning decal and release trigger to lock. The cow horn has 3 alternative positions for lifting and a further one for the storage/fold position.

### Cow Horn Lift Positions





### Lifting from a sitting position

Select the most suitable cow horn position to suit the patient or environment. Once complete, position the Standaidd so that the patient's feet fit in the foot tray and the knees rest against the kneepad. The optimal position for the kneepad is just below the patella. Adjust the kneepad accordingly and fit the safety belt if required. If not required, the safety belt can be detached and stored away.

Attach the loops of the sling to the grey hooks at the end of the Oxford Journey handgrips. The patient's hands can be placed on the handgrips (see Fig A or C) or they may also grasp the sling straps or trailing loops if it's more comfortable to do so, avoiding any potential pinch points. Prior to lifting make a visual check to ensure the loops of the sling are secure and in place. Then, push the UP button on hand control to raise the patient to the standing position (see Fig B or D).


**Standing Sling:****Fig. A****Fig. B****Transport Sling:****Fig. C****Fig. D**

**NOTE:** Oxford recommends a risk assessment be completed before carrying out your lifting operation.

**⚠ WARNING**

If using the Oxford Journey for gait training, the foot tray should be removed and a thorough risk assessment carried out to ascertain suitability prior to lifting.

## 6. Smart Monitor & Handset

## Intended Use

The Oxford Smart Monitor is a control system for Oxford Hoists. The Smart Monitor stores useful servicing information about the hoist that can be recalled when required. This servicing information includes:

- Number of patient lift cycles
- Total work done by actuator
- Number of lift overloads (attempted lifts above the safe working load)
- Number of days since last service interval

The Oxford Smart Monitor contains a microprocessor inside making it possible to read out service data via the on-board LCD screen.

## Features

- Data collection in the control box with data display via LCD screen
- Work counter
- Intelligent cycle counter
- Service indicator
- Service interval indicator
- Overload information
- 3 step battery indicator

## Benefits

- Improved safety for both patient and carer
- Accurate service data available at the touch of a button
- Optimised product life time
- Ease of maintenance for engineers and service technicians
- Fast and easy error identification

## LCD Display Screen

The Oxford Smart Monitor has the option to read out information via the LCD display screen. It is possible to read out total lifting cycles, total work done, overloads and number of days since last service, which can be used to quickly and easily evaluate the condition of the lift actuator.

This information is accessed by a ½ second press on the 'UP' button on the User handset or redundant 'UP' button on the Smart Monitor.

|                                                                                   |         |                                                                            |
|-----------------------------------------------------------------------------------|---------|----------------------------------------------------------------------------|
|  | 12034   | Total Lifting Cycles                                                       |
|  | 1257000 | Total Work Done (Amps * Seconds)                                           |
|  | 7       | Number of Overloads                                                        |
|  | 90/180  | Number of Days (since last service) /<br>Number of Days (between services) |

## Service & Usage Information

When it is time for a service, the service symbol will appear on the display. In addition, when the user handset is activated, the Smart Monitor will give an audible signal giving notice to users that a service is required.

The service symbol will be displayed each time the raise or lower keys are depressed on the User handset. When the service symbol is displayed, contact your authorised service provider immediately to arrange a service.


If the lift stops because of an overload (an attempt to lift more than the safe working load) the overload symbol will appear on the LCD display and the lift will cease to operate until the additional load (above the SWL) has been removed.

### **WARNING**

**If the LCD display screen on your Smart Monitor has recorded an overload, Oxford recommend that ALL routine daily checks are carried out on the lift prior to further use (for a list of daily checks, please refer to your lift user manual). Oxford also recommend that you contact your authorised Oxford service provider for additional guidance.**

## On-Board Redundant Controls


On-board redundant controls enable the lift to be raised or lowered in the event of an emergency.


## Battery Information

The display showing full battery means that the battery is fully charged and the lift is ready for use.

The display showing a half empty battery indicates that it is time to charge the battery.

The empty battery symbol showing on the display indicates that the battery has no capacity left and should be placed on charge immediately.

### **WARNING**

To avoid possible permanent damage to the battery, the battery should be placed on charge as soon as the display indicates the half empty battery symbol.

## 7. Operating Instructions

### **WARNING**

**Oxford patient lifts should only be operated by fully trained/competent care givers.**

#### 1. Leg Adjustment

The legs on the Oxford Journey are adjustable for width. The legs can be opened to enable access around armchairs or wheelchairs. For transferring and negotiating narrow doorways/passages, the lift legs should be in the closed position. To achieve the adjustment, the leg adjuster pedal, located at the rear of the base, is compressed right (DOWN) to open the legs outwards and left (UP) to close the legs. The adjustment can be carried out with the patient in the lift, but whether loaded or unloaded the adjustment should be made when the lift is moving.

#### 2. Castors and Braking

The lift has two braked castors, which can be applied for parking. When lifting, the castors should be left free and un-braked unless an appropriate and thorough risk-assessment determines otherwise. The lift will then be able to move to its natural centre of gravity. If the brakes are applied, it is the patient that will swing to the centre of gravity and this may prove disconcerting and uncomfortable.

An example of when brakes might be applied is if a patient is likely to push the lift away with their feet prior to lifting, risking injury to patient and/or carer. Another example is if the stand-aid is being used without the footplate (patient's feet on the floor), in this case, applying the brakes will ensure the lift doesn't 'run away' from the patient when being raised, leading to an avoidable fall and possible injury to patient and/or carer.

#### 3. Raising & Lowering the Boom

A powerful electric actuator that is controlled by a simple hand control unit achieves the movement of the boom. The hand control has two buttons with directional arrows UP and DOWN ↑↓. The actuator stops automatically at its limit of travel in both directions. The hand control plugs into a socket at the base of the control box.

#### 4. Emergency Stop

The red Emergency Stop Button is located on the front of the control box and is activated by pressing it in. This will cut all power to the lift and will only be reset by twisting the button clockwise and releasing.

#### 5. Redundant Controls

All Oxford lifts are fitted with raise and lower buttons on the control box. These are located underneath the emergency stop button and can be used to lower/raise the patient should the hand control fail.

#### 6. Mechanical Emergency Down

In the case of a complete electrical failure the electrical actuator is fitted with mechanical lowering device (RED BOSS). This will only operate when the lift is under load. The device must be pulled upwards to activate, and a slow descent will commence.

**Note:** The MINIMUM load required to manually lower the hoist using the mechanical emergency down function is 30 kg.

## 7. Batteries

A LOW VOLTAGE ALARM protects the batteries from deep discharge. This will sound when the batteries need recharging and the hand control is being operated. It will not sound independently of the hand control being operated. **DO NOT IGNORE THIS WARNING ALARM.** Complete the lifting operation and place the battery on charge (see charging instructions). The battery charge level is displayed by the indicator on the control box.

## ⚠ CAUTION

**The battery, charger, hand control, control box and actuator should not be opened by unauthorised personnel. (Contact your distributor for warranty and repairs).**


## 8. Slings

The Oxford Journey is suitable for patients in the SITTING, position. The slings suitable for this device are listed as follows:

- Oxford Transport sling
- Oxford Standing sling

## 9. Knee Pad Adjustment


The Oxford Journey has an adjustable kneepad. It is recommended for safety and ergonomic reasons that any adjustment is carried out by holding the top of the adjustment handle. The range of travel is 10.2 inches (260mm). Adjustment is achieved by following two easy steps. STEP 1 turn the kneepad adjustment lever to the left to unlock. STEP 2 either push handle downwards or pull handle upwards to select correct position. Once the ideal position is found, simply turn the kneepad adjustment lever to the right to lock in place, NOTE: the lever can be operated from either side of the lift.

## 10. Removable footplate

The Oxford Journey has a removable footplate. To remove the footplate, move to the front of the Oxford Journey and lift the footplate upward from the middle of the base assembly.

**NOTE:** Be careful when removing the footplate from its low-seated position as it is heavy. Always use appropriate moving and handling techniques.


## **⚠ WARNING**


If using the Oxford Journey for gait training, the foot tray should be removed and a thorough risk assessment carried out to ascertain suitability prior to lifting.

## 11. Adjustable cow horn


The Oxford Journey has an adjustable cow horn. This allows the carer to position the cow horn in a comfortable position to suit the patient or environment. The cow horn has 3 positions. Position 1 is more suitable for a large/taller patient, position 2 is more suitable for a medium/average size patient and position 3 is more suitable for a shorter/pediatric patient or when standing a patient from a low chair. (Please note, these are guidelines only, it may be necessary to adjust according to patient comfort). To adjust the cow horn, simply pull the trigger backwards towards the rear of the lift whilst holding the cow horn. Then select the desired position as indicated by the positioning decal and release trigger to lock. The cow horn has 3 alternative positions for lifting and a further one for the storage/fold position.


**Position 1**

**Position 2**

**Position 3**

## **Cow Horn Lift Positions**

## 8. Safety Precautions

Please read and follow the safety precautions listed below. The operation and use of Oxford patient lifts is simple and straightforward. Following these few basic safety precautions will make lifting operations easy and trouble free.

### **READ AND UNDERSTAND THE USER INSTRUCTION MANUAL BEFORE USING YOUR “JOURNEY” LIFT**

**WARNING:** Important safety information for hazards that might cause serious injury.

**CAUTION:** Information for preventing damage to the product.

**NOTE:** Information to which you should pay special attention.

### **WARNING**

- Oxford patient lifts should only be operated by full trained/competent care givers.
- **ALWAYS** plan your lifting operations before commencing.
- **ALWAYS** carry out the **DAILY CHECK LIST** before using the lift.
- **ALWAYS** familiarise yourself with the operating control and safety features of a lift before lifting a patient.
- **DO NOT** use a sling unless it is recommended for use with the lift.
- **ALWAYS** check the sling is suitable for the particular patient and is of the correct size and capacity.
- **NEVER** use a sling, which is frayed or shows signs of damaged.
- **ALWAYS** fit the sling according to the instructions provided (user instructions).
- **DO NOT** exceed the safe working load.
- **ALWAYS** check the safe working load of the lift is suitable for the weight of the patient.
- **ALWAYS** carry out lifting operations according to the instructions in the user manual.
- **NEVER** disconnect or bypass a control or safety feature because it seems easier to operate the lift.
- **DO NOT** attempt to manoeuvre the lift by pushing on the mast, boom or patient.
- **ALWAYS** manoeuvre the lift with the handle / foot push pad provided.
- **ALWAYS** lower the patient to the lowest comfortable position before transfers.
- **DO NOT** push a loaded lift at speed, which exceeds a slow walking pace (3 Km/hour 0.8 metres/second).
- **DO NOT** push the lift over uneven or rough ground. Particularly if loaded.
- **DO NOT** attempt to push/pull a loaded lift over a floor obstruction.
- **NEVER** force an operating/safety control. All controls are easy to use and do not require excessive force.
- The **MINIMUM** load required to manually lower the hoist using the mechanical emergency down function is 30 kg.

- **DO NOT** park a loaded lift on ANY sloping surface.
- **DO NOT** use electric lifts in a shower.
- **DO NOT** charge an electric lift in a bathroom or shower room.
- **DO NOT** lift a patient unless you are trained and competent to do so.
- **DO NOT** use the lift, or allow the lift to be used, for any other purpose.
- **DO NOT** bump the lift down steps, loaded or unloaded.
- **DO NOT** attempt to negotiate a loaded lift on a slope, which exceeds 1:12 (approximately 5 degrees).
- **DO NOT** attempt to negotiate a slope without a second helper being present.
- **DO NOT** use in a wet or corrosive environment such as poolside locations.

## 9. Charging Instructions

When the power pack needs charging it is removed from the lift and fitted to an off board charging unit.

Oxford recommend an additional battery pack is purchased, so that one pack can be on charge at all times.

1. Remove the power pack from the lift. The pack is retained by a simple latch at the top of the power pack. Lift the latch and the power pack will be released.
2. Fit the power pack to the charging unit. The location and latching of the power pack to the charger is the same system as used on the lift.
3. Plug the charger mains plug into a suitable mains outlet and switch the mains supply ON.
4. Charging is automatic and will fully charge the batteries over a period of four (4) hours.  
**Note:** Even if the charger is left plugged in for extended periods it will not allow the batteries to overcharge.
  - a) Green Light - Indicates main power is on.
  - b) Yellow Light - Indicates battery is charging.
  - c) Battery will be fully charged when yellow light goes off.
5. To return the lift to service, switch OFF the mains supply and remove the power pack from the charger. Fit the power pack to the lift and make sure the latch holding the pack in place is fully engaged. "Click" in place. The charging of Oxford electric lifts is simple and straightforward, but it is important to follow the charging instructions closely. Please pay particular attention to the following points, they will help you avoid problems with discharged batteries.

### WARNING

- **The battery, charger, hand control & control box are NOT to be opened by unauthorised personnel. (Contact your distributor for warranty and repairs).**
- **DO NOT touch battery/charger terminals.**
- **KEEP the batteries fully charged. Place the battery on charge whenever it is not in use. If it is more convenient to do so, place on charge every night. The charger will not allow the batteries to overcharge.**
- **NEVER run the batteries completely flat. As soon as the audible warning sounds, complete the lifting operation in hand and place on charge.**
- **To avoid possible permanent damage to the battery, the battery should be placed on charge as soon as the display indicates the half empty battery symbol.**
- **NEVER store the battery for long periods without regular charging throughout the storage period.**
- **ALWAYS make sure the mains power to the charger is switched off before connecting or disconnecting the battery pack.**
- **NEVER leave the battery pack connected with the charger mains power switched off.**

- **DO NOT** leave the charger switched on with the battery disconnected.
- **DO NOT** use the lift while charging is taking place.
- **CHECK** the lift is not charging before moving as the electrical connection may be damaged.
- **NEVER** disconnect the charger plug by pulling on the cable.
- **BE CAREFUL** not to trip over the charge lead.
- **DO NOT** charge an electric lift in a bathroom or shower room.

## 10. Maintenance Schedule & Daily Check List

All Oxford products are designed for minimum maintenance, however some safety checks and procedures are required. Schedules of DAILY tasks are detailed below. Daily checks and a six monthly service, inspection and test will ensure a lift is kept in optimum safe working condition. A list of spare parts is available upon request.

The **LOAD TEST** and certification should only be carried out by qualified personnel or an authorised service agent / dealer.

### DAILY CHECK LIST:

Oxford strongly recommends the following checks be carried out on a daily basis and before using lift.

- ENSURE the lift moves freely on its castors.
- EXAMINE the spreader bar (cow horns) is firmly attached to the boom.
- EXAMINE the sling hooks/clips on the spreader bar for damage and/or excessive wear.
- ENSURE the legs open and close correctly and lock in all open and closed positions.
- OPERATE the hand control to confirm the boom raises and lowers satisfactorily.
- ENSURE the control box display panel DOES NOT have a spanner symbol illuminated and no audible beep can be heard when depressing raise/lower buttons.
- CHECK the operation of the emergency stop button.
- EXAMINE slings for fraying or other damage. DO NOT use any sling if damaged or if the sling shows signs of wear.
- CONFIRM the lift is not sounding a low battery alarm when the hand control is operated. If the alarm sounds, DO NOT use, and place on charge immediately.
- MAKE SURE the actuator connection points are secure.
- VISUALLY INSPECT all fixings and fasteners and ensure they are secure prior to use.

### MAINTENANCE, INSPECTION AND TEST

Oxford recommends a thorough inspection and test of the Oxford Journey lift and lifting accessories; slings, etc. are carried out on a regular basis. Inspection frequency varies depending on locality, so you must check with your dealer or local government agency as appropriate regarding how often an inspection is required. The examination and test should be conducted according to the recommendations and procedures provided in this manual. Whenever possible Oxford recommends maintenance, inspection and certified testing is carried out by authorised service agent / dealers only.

**NOTE:** These recommendations are in compliance with the requirements of 1998 No2307 Health and Safety: The Lifting Operations and Lifting Equipment Regulations 1998. This is a UK regulation. Outside the UK please check your local country requirements.

### Cleaning

Clean with ordinary soap and water and/or any hard surface disinfectant. Harsh chemical cleaners or abrasives should be avoided as these may damage the surface finish of the lift. Avoid wetting any of the electrical parts. After cleaning, the unit should be thoroughly dried.

### End of Life Disposal

**Battery:** All batteries in the product and any accompanying accessories must be recycled separately. Batteries should be disposed of in accordance with local or national regulations.

**Slings:** Slings and associated material accessories should be sorted as combustible waste in accordance with local or national regulations.

**\*THESE CHECKS SHOULD INCLUDE:**

|                                                                                                                                                                                                                                                                                                                                             | Initially | Before Use | Service Intervals |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|------------|-------------------|
| <b>1. BOOM:</b> Check the attachment of the boom to the mast. Make sure there is only minimal side movement of the boom and the boom is free to rotate on the boom bearing.                                                                                                                                                                 | ✓         | ✓          | ✓                 |
| <b>2. SLING ATTACHMENT POINTS:</b> Examine the sling hooks and attachment clips on the spreader bar and boom for excessive wear.                                                                                                                                                                                                            | ✓         | ✓          | ✓                 |
| <b>3. MAST:</b> Check the operation of the mast-locking device. Make sure the mast fully engages into the socket.                                                                                                                                                                                                                           | ✓         | ✓          | ✓                 |
| <b>4. COW HORN:</b> Check cow horn is engaged and in correct position.                                                                                                                                                                                                                                                                      | ✓         | ✓          | ✓                 |
| <b>5. ACTUATOR:</b> The actuator should require no maintenance other than checking for correct operation and listening for unusual noise.                                                                                                                                                                                                   |           |            | ✓                 |
| <b>6. KNEE PAD:</b> Check knee pad adjuster is fastened securely.                                                                                                                                                                                                                                                                           | ✓         | ✓          |                   |
| <b>7. SMART MONITOR:</b> Check the function of the emergency stop button. Inspect the hand control socket for correct fitting. Check functioning of the hand control. Check the redundant controls and confirm they operate as intended. The MINIMUM load required to manually lower the hoist using the emergency lower function is 30 kg. | ✓         | ✓          | ✓                 |
| <b>8. BATTERIES:</b> The batteries are housed in the power pack and should not require maintenance other than the regular charging as detailed in the charging instructions. Check that the connections remain clean.                                                                                                                       |           |            | ✓                 |
| <b>9. BATTERIES (CHARGING):</b> Confirm the lift is not giving a low battery alarm when the hand control is operated. If the alarm sounds, do not use and place the battery on charge immediately.                                                                                                                                          | ✓         | ✓          |                   |
| <b>10. LEG ADJUSTMENT:</b> Check the legs operate in both full extensions (inward/outward).                                                                                                                                                                                                                                                 | ✓         | ✓          | ✓                 |
| <b>11. CASTORS:</b> Check all castors for firm attachment to the legs. Check for free rotation of the castor and the wheels.                                                                                                                                                                                                                | ✓         | ✓          | ✓                 |
| <b>12. CLEANING:</b> Using a damp cloth, wipe clean with ordinary soap and water and/or any hard surface disinfectant. Harsh chemical cleaners or abrasives should be avoided as these may damage the surface finish of the lift. Avoid wetting any of the electrical parts. After cleaning, the unit should be thoroughly dried.           |           | ✓          | ✓                 |
| <b>13. LOAD TEST:</b> The load test should be carried out in accordance with the manufacturers test procedures. It is strongly recommended that an authorised service dealer carry out the test.                                                                                                                                            |           |            | ✓                 |
| <b>14. BASE AND WHEELS:</b> Ensure base is even and level (all four wheels are on the floor).                                                                                                                                                                                                                                               | ✓         | ✓          | ✓                 |
| <b>15. SLINGS:</b> Check for wear and fraying.                                                                                                                                                                                                                                                                                              | ✓         | ✓          | ✓                 |
| <b>16. LUBRICATION:</b> Lubricate pivot joints with a silicone based oil, including mast and boom connections, pedal assembly, spreader bar joint (only if required).                                                                                                                                                                       |           |            | ✓                 |
| <b>17. HAND SET:</b> Ensure plugged fully into controller.                                                                                                                                                                                                                                                                                  | ✓         | ✓          |                   |

# 11. Technical Specifications

|                                              |                   |         |
|----------------------------------------------|-------------------|---------|
| Safe Working Load .....                      | 341 lbs.....      | 155 kgs |
| Maximum Overall Length .....                 | 39.4 inches.....  | 1000 mm |
| Minimum Overall Length .....                 | 36.8 inches ..... | 935 mm  |
| Maximum Overall Height.....                  | 63 inches .....   | 1600 mm |
| Minimum Overall Height.....                  | 44.5 inches ..... | 1130 mm |
| Maximum Height to attachment point.....      | 60.2 inches ..... | 1530 mm |
| Minimum Height to attachment point.....      | 28.4 inches ..... | 720 mm  |
|                                              |                   |         |
| Turning Radius .....                         | 45.1 inches.....  | 1145 mm |
| Legs Open - External Width .....             | 38.2inches .....  | 970 mm  |
| Legs Open - Internal Width .....             | 34.5 inches ..... | 875 mm  |
| Legs Closed - External Width .....           | 27.0 inches ..... | 685 mm  |
| Legs Closed - Internal Width .....           | 22.2 inches ..... | 590 mm  |
| Widest Point (between support handles) ..... | 26.8 inches ..... | 680 mm  |
|                                              |                   |         |
| Overall Height of Legs.....                  | 4.3 inches .....  | 110 mm  |
| Ground Clearance .....                       | 1.2 inches .....  | 30 mm   |
| Front Twin Castors .....                     | 3.0 inches .....  | 75 mm   |
| Rear Braked Castors.....                     | 4.0 inches .....  | 100 mm  |

## Weights

|                                      |               |         |
|--------------------------------------|---------------|---------|
| Mast, Base & Boom Assembly .....     | 82.4 lbs..... | 37.4 kg |
| Power Pack .....                     | 6.6 lbs.....  | 3.0 kg  |
| Total .....                          | 89.0 lbs..... | 40.4 kg |
| Base Assembly (not inc battery)..... | 30.2 lbs..... | 13.7 kg |
| Mast & Boom (not inc battery) .....  | 42.3 lbs..... | 19.2 kg |
| Foot Tray .....                      | 9.9 lbs.....  | 4.5 kg  |

**NOTE:** Manufacturing tolerances apply to all dimensions.

## **ELECTRICAL SPECIFICATIONS:**

|                            |                                            |
|----------------------------|--------------------------------------------|
| BATTERIES .....            | 24 volt Rechargeable sealed lead acid type |
| BATTERY CAPACITY .....     | 3.2 Ampere hours                           |
| CHARGER RATED INPUT .....  | 100 - 240 V AC/ 50/60 Hz                   |
| CHARGER RATED OUTPUT ..... | 29.5 VDC, Max. 19 W                        |

## **ELECTRIC SHOCK PROTECTION**

|               |                       |
|---------------|-----------------------|
| CHARGER ..... | CLASS II (EN 60601-1) |
| LIFT .....    | INTERNAL POWER SOURCE |

## **DEGREE OF SHOCK PROTECTION**

|               |                     |
|---------------|---------------------|
| CHARGER ..... | TYPE B (EN 60601-1) |
| LIFT .....    | TYPE B (EN 60601-1) |

## **ENVIRONMENTAL CONDITIONS:**

Outside this environment functionality and safety may be compromised.

### **OPERATING:**

|                            |                                    |
|----------------------------|------------------------------------|
| TEMPERATURE .....          | 5°C to 40°C                        |
| RELATIVE HUMIDITY .....    | 20% to 90% @ 30°C - not condensing |
| ATMOSPHERIC PRESSURE ..... | 700 to 1060 hPa                    |

### **STORAGE:**

|                            |                                    |
|----------------------------|------------------------------------|
| TEMPERATURE .....          | -10°C to +50°C                     |
| RELATIVE HUMIDITY .....    | 20% to 90% @ 30°C - not condensing |
| ATMOSPHERIC PRESSURE ..... | 700 to 1060 hPa                    |

### **IP RATINGS**

|                         |      |
|-------------------------|------|
| CONTROL BOX .....       | IPX4 |
| ACTUATOR .....          | IPX4 |
| OFF BOARD CHARGER ..... | IPX5 |
| BATTERY .....           | IPX5 |
| HANDSET .....           | IPX5 |

### **DUTY CYCLES**

|                |                      |
|----------------|----------------------|
| ACTUATOR ..... | 10% (2 min./18 min.) |
| CHARGER .....  | Approx. 4 hrs.       |

### **ACOUSTICS**

|                                    |           |
|------------------------------------|-----------|
| A-Weighted sound power level ..... | 60 dB (A) |
|------------------------------------|-----------|

**KEY SYMBOLS:**

The following symbols are used on the charger, control unit and battery:

This product is a Medical Device in accordance with EU Medical Device Regulation 2017/745.

Type B equipment, as per EN 60601-1.

Class 2 equipment.

The disposal of batteries should not be mixed with general household waste. Batteries should be disposed of in accordance with local or national regulations.

The disposal of electronics should not be mixed with general household waste.

For indoor use.

ATTENTION, consult accompanying documents.

Refer to the user instructions before use (Blue background).

Certification mark indicating conformity with the applicable requirements for products sold within Great Britain.

A point of contact between the non-European medical device manufacturer, the national Competent Authorities (Ministry of Health) and Notified Bodies.

A unique device identification intended to provide single, globally harmonised positive identification of medical devices through distribution and use.

Indicates the entity importing the medical device into the locale.

## 12. Warranty

Oxford has an established network of authorised distributors and service agents who will be pleased to handle all of your purchasing, warranty, repair and maintenance enquiries. It is recommended that our products are commissioned by your authorised distributor or service agent.

The distributor or service agent operates the warranty programme, so it is important to keep a record of their name, address and telephone number so they can be contacted should any problem arise. If you are in any doubt as to where your lift was purchased, Oxford can trace the supplier if you quote the serial number of the lift.

All warranty claims are subject to the correct use and maintenance in accordance with the user instruction manual provided with each lift. Damage caused by use in unsuitable environmental conditions or failure to maintain the product in accordance with user and service instructions is not covered. Your statutory rights remain unaffected.

**REMEMBER:** Contact your authorised distributor or service agent for purchases, repairs, servicing and certified maintenance.

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Oxford Medical Healthcare Ltd