

Data management/Device compatibility list

English

This supplement specifies the devices compatible with ResMed data management products. For full details regarding the correct use of these devices, see the relevant User Guides.

Note: Not all products are available in all regions.

	ResScan™ 7.1	ResScan™ Essentials 1.1	EasyCare Tx 9.0/ Tx Link™	EasyCare Tx 2 1.3/ Tx Link™ 2	AirView™	RCM	RCMH ⁱ
AirMini™					✓		
AirSense 11/AirCurve 11 Series							
AirSense™ 11 AutoSet	✓	✓		✓	✓		
AirSense™ 11 Elite	✓	✓		✓	✓		
AirSense™ 11 CPAP	✓	✓		✓	✓		
AirCurve™ 11 ASV PaceWave™	✓	✓		✓	✓		
AirCurve™ 11 S	✓	✓		✓	✓		
AirCurve™ 11 ST	✓	✓		✓	✓		
AirCurve™ 11 VAuto	✓	✓		✓	✓		
AirSense 10/AirCurve 10 Series							
AirSense™ 10 AutoSet/ AirSense™ 10 AutoSet for Her	✓	✓	✓		✓		
AirSense™ 10 Elite	✓	✓	✓		✓		
AirCurve™ 10 CS PaceWave™	✓	✓	✓		✓		
AirCurve™ 10 CS-A PaceWave™	✓	✓	✓		✓		
AirCurve™ 10 S	✓	✓	✓		✓		
AirCurve™ 10 ST	✓	✓	✓		✓		
AirCurve™ 10 VAuto	✓	✓	✓		✓		
Lumis Series							
Lumis™ 100 VPAP S	✓	✓	✓		✓		
Lumis™ 100 VPAP ST	✓	✓	✓		✓		
Lumis™ 150 VPAP ST	✓	✓	✓		✓		
Lumis™ 100 VPAP ST-A	✓	✓	✓		✓		
Lumis™ 150 VPAP ST-A	✓	✓	✓		✓		
Lumis™ HFT					✓		
Lumis™ ST EU 30cmH2O PAC					✓		
S9 Series							
Auto 25 (S9)	✓		✓				
AutoSet CS-A (S9 PaceWave™)	✓		✓				
AutoSet CS™ (S9 PaceWave™) ⁱⁱ	✓		✓		✓		
AutoSet CS™ (S9) ⁱⁱ	✓		✓		✓		
S9 AutoSet™ ⁱⁱ	✓		✓		✓		
S9 Elite™ ⁱⁱ	✓		✓		✓		
S9 Escape™	✓		✓				
S9 Escape™ Auto	✓		✓				
S9 VPAP™ Tx			✓				
V-Auto (S9) ⁱⁱ	✓		✓		✓		
VPAP™ Adapt (S9 PaceWave™)	✓		✓				
VPAP™ Adapt (S9)	✓		✓				
VPAP™ S (S9)	✓		✓				
VPAP™ ST (S9) ⁱⁱ	✓		✓		✓		



	ResScan™ 7.1	ResScan™ Essentials 1.1	EasyCare Tx 9.0/ Tx Link™	EasyCare Tx 2 1.3/ Tx Link™ 2	AirView™	RCM	RCMH ⁱ
VPAP™ ST (S9 iVAPS) ⁱⁱ	✓		✓		✓		
VPAP™ ST-A (S9 iVAPS)	✓		✓				
VPAP™ ST-A (S9) ⁱⁱ	✓		✓		✓		
S8 Series							
S8 AutoScore™	✓		✓				
S8 AutoScore™ II	✓		✓				
S8 AutoSet Spirit™	✓		✓				
S8 AutoSet Spirit™ II	✓		✓				
S8 AutoSet Vantage™	✓		✓				
S8 AutoSet™ II	✓		✓				
S8 Elite™	✓						
S8 Elite™ II	✓		✓				
S8™ Auto 25	✓		✓				
VPAP™ IV Series	✓		✓				
S7 Series							
AutoSet CS™ 2	✓		✓				
S7 AutoSet Spirit™	✓						
S7 Elite™	✓						
VPAP™ Adapt SV	✓						
VPAP™ II Series	✓						
VPAP™ III Series	✓		✓				
VPAP™ III ST-A with QuickNav	✓		✓				
Ventilator							
Astral™ 100	✓	✓			✓	✓	✓
Astral™ 150	✓	✓			✓	✓	✓
Stellar™ 100	✓	✓	✓		✓	✓	
Stellar™ 130	✓	✓	✓		✓	✓	
Stellar™ 150	✓	✓	✓		✓	✓	
Elisée™ 150	✓						
VS III™	✓						
Home sleep testing device							
ApneaLink™ Air					✓		

RCM — ResMed Connectivity Module
RCMH — ResMed Connectivity Module Hospital

ⁱFor a list of available data for use by EMR system, please contact your ResMed representative.

ⁱⁱCompatible with AirView via SD card uploads only.



Made in USA
Manufacturer: ResMed Corp 9001 Spectrum Center
Boulevard San Diego CA 92123 USA



Made in Australia
Manufacturer: ResMed Pty Ltd 1 Elizabeth Macarthur
Drive Bella Vista NSW 2153 Australia

See ResMed.com for other ResMed locations worldwide. AirMini, AirSense, AirCurve, AirView, Astral, ApneaLink, Tx Link, TxLink 2, Lumis, S8, S9, Stellar, VPAP, AutoSet, S9 Escape, S9 Elite, S9 Escape Auto, S8 Elite, S8 AutoSet Vantage, PaceWave, ResScan, VPAP Tx, S8 AutoSet, S9 AutoSet, S9 Auto 25, VS III, Elisée, VPAP Adapt SV, AutoSet CS, S8 AutoScore, S8 AutoSet Spirit and VPAP are trademarks and/or registered trademarks of the ResMed family of companies. For patent and other intellectual property information, see ResMed.com/ip.

1018030/29 2025-08