

A photograph of a patient lying in a hospital bed, wearing a Philips AF541 oro-nasal mask. The mask is secured with a black headgear. A healthcare professional wearing blue gloves is adjusting the mask. The patient is looking towards the camera. The background is a soft-focus hospital room.

PHILIPS

AF541

Oro-nasal mask

Reduce skin breakdown while providing patient comfort

Philips AF541 oro-nasal mask for single-patient use

Delivering high-quality respiratory care to your patients in today's complex healthcare environment demands new thinking — especially when it comes to patient comfort.

Designed to rest comfortably on the face, the Philips AF541 non-invasive ventilation (NIV) mask offers interchangeable under-the-nose (UTN) and over-the-nose (OTN) cushions to help reduce skin breakdown at mask contact points.

Key advantages

- Available for patients > 40 lb (20 kg)
- Available in both CapStrap and four-point headgear
- All AF541 click-style mask elbows include nebulizing port
- Can be used with interchangeable elbows
- CleanClip system provides convenient storage for mask when not in use

Ordering information



AF541 with over-the-nose cushion (OTN) can be easily lifted off the face or rotated to free up hands during oral care.

AF541 with over-the-nose cushion (OTN), single-patient use

AF541 with CapStrap Headgear	Size	Quantity	Product #
AF541 with EE Leak 1 elbow	Small	10	73-1120915
AF541 with EE Leak 1 elbow	Medium	10	73-1120916
AF541 with EE Leak 1 elbow	Large	10	73-1120917
AF541 with EE Leak 1 elbow	Extra Large	10	73-1120918
AF541 with EE Leak 2 elbow	Small	10	73-1120935
AF541 with EE Leak 2 elbow	Medium	10	73-1120936
AF541 with EE Leak 2 elbow	Large	10	73-1120937
AF541 with EE Leak 2 elbow	Extra Large	10	73-1120938
AF541 with SE Standard elbow	Small	10	73-1120945
AF541 with SE Standard elbow	Medium	10	73-1120946
AF541 with SE Standard elbow	Large	10	73-1120947
AF541 with SE Standard elbow	Extra Large	10	73-1120948
AF541 with Four-Point Headgear	Size	Quantity	Product #
AF541 with EE Leak 1 elbow	Small	10	73-1120910
AF541 with EE Leak 1 elbow	Medium	10	73-1120911
AF541 with EE Leak 1 elbow	Large	10	73-1120912
AF541 with EE Leak 1 elbow	Extra Large	10	73-1120913
AF541 with EE Leak 2 elbow	Small	10	73-1120930
AF541 with EE Leak 2 elbow	Medium	10	73-1120931
AF541 with EE Leak 2 elbow	Large	10	73-1120932
AF541 with EE Leak 2 elbow	Extra Large	10	73-1120933
AF541 with SE Standard elbow	Small	10	73-1120940
AF541 with SE Standard elbow	Medium	10	73-1120941
AF541 with SE Standard elbow	Large	10	73-1120942
AF541 with SE Standard elbow	Extra Large	10	73-1120943

Not all products are available in all markets.

Ordering information *(continued)*



AF541 with under-the-nose cushion (UTN) can be used to offload mask pressure at bridge of nose.

AF541 with under-the-nose cushion (UTN), single-patient use

AF541 with CapStrap Headgear	Size	Quantity	Product #
AF541 with EE Leak 1 elbow	Extra Small	10	73-1133660
AF541 with EE Leak 1 elbow	A	10	73-1133661
AF541 with EE Leak 1 elbow	B	10	73-1133662
AF541 with EE Leak 1 elbow	C	10	73-1133663
AF541 with EE Leak 2 elbow	Extra Small	10	73-1133676
AF541 with EE Leak 2 elbow	A	10	73-1133677
AF541 with EE Leak 2 elbow	B	10	73-1133678
AF541 with EE Leak 2 elbow	C	10	73-1133679
AF541 with SE Standard elbow	Extra Small	10	73-1133684
AF541 with SE Standard elbow	A	10	73-1133685
AF541 with SE Standard elbow	B	10	73-1133686
AF541 with SE Standard elbow	C	10	73-1133687
AF541 with Four-Point Headgear	Size	Quantity	Product #
AF541 with EE Leak 1 elbow	Extra Small	10	73-1133626
AF541 with EE Leak 1 elbow	A	10	73-1133627
AF541 with EE Leak 1 elbow	B	10	73-1133658
AF541 with EE Leak 1 elbow	C	10	73-1133659
AF541 with EE Leak 2 elbow	Extra Small	10	73-1133672
AF541 with EE Leak 2 elbow	A	10	73-1133673
AF541 with EE Leak 2 elbow	B	10	73-1133674
AF541 with EE Leak 2 elbow	C	10	73-1133675
AF541 with SE Standard elbow	Extra Small	10	73-1133680
AF541 with SE Standard elbow	A	10	73-1133681
AF541 with SE Standard elbow	B	10	73-1133682
AF541 with SE Standard elbow	C	10	73-1133683

Not all products are available in all markets.

Ordering information *(continued)*



Leak 1 entrainment elbow (EE) click-style has an anti-asphyxia valve for dedicated non-invasive ventilators (use with single-limb circuits)



Leak 2 entrainment elbow (EE) click-style has an anti-asphyxia valve for dedicated non-invasive ventilators and intentional leak ports (use with single-limb circuits)



Standard elbow (SE) click-style is typically used on ventilators equipped with integrated safety valve systems (use with dual-limb circuits)



Bronchoscopy elbow click-style is compatible with Philips Respironics masks that use the click style, 22 mm elbow design

AF541 accessories and replacement parts

Elbows		Quantity	Product #
EE Leak 1 elbow with NIVO port		5	73-1120967
EE Leak 2 elbow with NIVO port		5	73-1120968
SE Standard elbow with NIVO port		5	73-1120969
Bronchoscopy elbow		5	73-1120984
Cushions, Over-The-Nose (OTN)	Size	Quantity	Product #
AF541 over-the-nose	Small	10	73-1120960
AF541 over-the-nose	Medium	10	73-1120961
AF541 over-the-nose	Large	10	73-1120962
AF541 over-the-nose	Extra Large	10	73-1120963
Cushions, Under-The-Nose (UTN)	Size	Quantity	Product #
AF541 under-the-nose	Extra Small	10	73-1127634
AF541 under-the-nose	A	10	73-1120980
AF541 under-the-nose	B	10	73-1120981
AF541 under-the-nose	C	10	73-1120982
Headgear	Size	Quantity	Product #
AF541 CapStrap	Small/Medium	10	73-1126500
AF541 CapStrap	Large/Extra Large	5	73-1120972
AF541 Four-Point with forehead pad	Small/Medium	5	73-1120970
AF541 Four-Point with forehead pad	Large/Extra Large	5	73-1120971

Not all products are available in all markets.

Tri-anim provides innovative respiratory, anesthesia, critical care and neonatal care products to hospitals, health systems and other patient care facilities nationwide.



Contact your local Tri-anim Account Manager
for more information or visit: tri-anim.com

Caution: U.S. federal law restricts these devices
to sale by or on the order of a physician.



800.874.2646 • www.tri-anim.com