

QA3 **DRIVE**TM

Emergency Stretcher



Specification Sheet

AneticAid[®] 

QA3 ///IDRIVE™

Emergency Stretcher

Inbuilt motorised drive

No need to push or pull the stretcher: our inbuilt motorised DRIVE feature linked to the stretcher's fifth wheel offers superb manoeuvrability and control

Pushing handles

Four ergonomically positioned pushing handles for supreme manoeuvrability

Effortless travel

The QA3™ DRIVE Emergency Stretcher can travel up to gradients of 6°.

X-ray translucent platform

Full length radiotranslucent platform - suitable for digital or traditional plates.

Long-life battery

Capable of **20 kilometers on a single charge**, the DRIVE feature is powered by an advanced nickel metal hydride rechargeable battery.

For more information on any of our products or service contracts, call:

sales@aneticaid.com **+44 (0)1943 878647** aneticaid.com



DESIGN & INNOVATION



MADE IN THE UK



GLOBAL DISTRIBUTION

Innovative **Medical** Technology - **Practically Applied**



Anetic Aid Ltd, Sapper Jordan Rossi Park, Otley Road, Baildon BD17 7AX, UK

Simultaneous braking low rolling resistance castors

Polyurethane tyres offer low rolling resistance for high manoeuvrability: 150mm diameter central total lock castors fitted as standard, helping to offer a rigid and reliable work platform.

Exceptional patient comfort

Low height of 535mm (excluding mattress) & store away side rails allow for easier mount/dismount, particularly for frail and elderly patients. Also including a versatile positioning headrest, wide patient platform and pressure care mattress for comfort.

Press on / Press off 5th wheel steer

Our unique mechanism ensures maximum manoeuvrability and means the QA3 is easily guided by one support worker.

Fixed IV pole

Innovative quick release safety locking mechanism, integrated into the head end side rail for easy access and storage.

Pressure care mattress

Manufactured from temperature-sensitive gel-feel foam that moulds to the patient's body shape, relieving pressure by distributing the load as evenly as possible. Outer material is anti-microbial, water resistant, vapour permeable and latex free. Now 5% resistant to chlorine and other disinfectant compounds.



SPECIFICATIONS		Code: 21122	
Minimum height	*535mm	Mattress length	2025mm
Maximum height	*835mm	Mattress width	700mm
*excludes mattress		Mattress depth	90mm
Minimum length	1820mm	Weight	157kg
Maximum length	2150mm	Backrest	0 to +90°
Overall width	805mm	Trendelenburg	0 to -12°
Castor diameter	150mm	Reverse Trendelenburg	*0 to -10°
Maximum patient weight		320kg	

*NOTE - Reverse Trendelenburg enabled through a build option

AneticAid®