

QA3™

Powered ⚡ Emergency Stretcher



Specification Sheet

AneticAid® 

QA3™

Powered ⚡ Emergency Stretcher

Powered functions

Four powered functions including platform raise and lower, back rest, Trendelenburg, Reverse Trendelenburg and CPR Return.

X-ray translucent platform

Full length radiotranslucent platform, suitable for digital or traditional plates.

Fixed transfusion pole

Innovative quick release safety locking mechanism, integrated into the head end side rail for easy access and storage.

Vertically operated side rails

Store away under the mattress, improving the ease of mounting and dismounting the stretcher for the patient.

Trendelenburg & Reverse Trendelenburg

Extremely smooth and controlled hydraulic mechanism, even at lowest height.

Simultaneous braking low rolling resistance castors

Polyurethane tyres offer low rolling resistance for high manoeuvrability: 150mm diameter central total lock castors fitted as standard, helping to offer a rigid and reliable work platform.

Exceptional patient comfort

Low height of 535mm (excluding mattress) & store away side rails allow for easier mount/dismount, particularly for frail and elderly patients. Also including a versatile positioning headrest, wide patient platform and pressure care mattress for comfort.

Press on / Press off 5th wheel steer

Our unique mechanism ensures maximum manoeuvrability and means the QA3 is easily guided by one support worker.

Pressure care mattress

Manufactured from temperature-sensitive gel-feel foam that moulds to the patient's body shape, relieving pressure by distributing the load as evenly as possible. Outer material is anti-microbial, water resistant, vapour permeable and latex free. Now 5% resistant to chlorine and other disinfectant compounds.



SPECIFICATIONS		Code: 21126	
Minimum height	*535mm	Mattress length	2025mm
Maximum height	*875mm	Mattress width	700mm
*excludes mattress		Mattress depth	90mm
Minimum length	1820mm	Maximum patient weight	320kg
Maximum length	2150mm	Backrest	0 to +70°
Overall width	805mm	KneeFlex	0 to +30°
Castor diameter	150mm	Trendelenburg	0 to -12°
Weight	162.5kg	Reverse Trendelenburg	0 to -10°

For more information on any of our products or service contracts, call:

sales@aneticaid.com **+44 (0)1943 878647** aneticaid.com



DESIGN & INNOVATION



MADE IN THE UK



GLOBAL DISTRIBUTION

Innovative **Medical** Technology - **Practically Applied**



Anetic Aid Ltd, Sapper Jordan Rossi Park, Otley Road, Baildon BD17 7AX, UK

AneticAid®